

Dialysate Instruments

PRICE LIST

APRIL 1, 2023

FDA Device Manufacturing License #s: 65387 (CA); 2020976 (Federal)



| MODEL | DIGITAL DIALYSATE METERS* | PRICE |
|--|--|-----------|
| D-4..... | Conductivity, Resistivity, TDS, Temperature..... | \$1984.00 |
| D-4BD..... | Conductivity, Resistivity, TDS, Temperature; bluDock™ enabled..... | \$2268.50 |
| D-6..... | Conductivity, Resistivity, TDS, pH, ORP, Temperature..... | \$2411.00 |
| D-6BD..... | Conductivity, Resistivity, TDS, pH, ORP, Temperature; bluDock enabled..... | \$2695.50 |
| *The D-6 and D-4 DIGITAL DIALYSATE METERS™ have received FDA 510(K) acceptance for use as medical devices for testing dialysate. | | |

| ITEM | ACCESSORIES | PRICE |
|--------------|--|----------|
| BLUDOCK..... | bluDock Wireless Data Transfer Module and instructions..... | \$284.50 |
| PKD..... | Hard foam-lined protective carry case with Standard/Buffer Solutions (large)*..... | \$265.00 |
| UPP..... | Hard foam-lined protective carry case (small)**..... | \$106.50 |
| UCC..... | Soft protective case with belt clip (Blue)..... | \$58.00 |
| UCCDT..... | Soft protective case with belt clip (Desert Tan)..... | \$58.00 |
| RPR..... | Replacement pH/ORP sensor..... | \$276.00 |
| DPDEKIT..... | Free Chlorine Measuring Kit..... | \$50.00 |

* Approximate dimensions: 13.6 L x 10 W x 5.4 D in./ 346 L x 254 W x 137 D mm

** Approximate dimensions: 9.5 L x 6.8 W x 3.6 D in./ 241 L x 171 W x 92 D mm

| MODEL | ANALOG DIALYSATE METERS | PRICE |
|----------|--|----------|
| D-1..... | Dialysate Meter™: 11-15 millimhos..... | \$497.00 |
| D-2..... | Dialysate Meter™: 0-10 & 0-20 millimhos..... | \$706.00 |

| ITEM | CALIBRATION STANDARDS - NIST Traceable | PRICE |
|---|---|----------------|
| Recommended Cond/TDS Standards are NaCl-14.0 and 442-3000. See Standard Solutions and Buffers Application Bulletin for additional values. Exact solution value must be specified when ordering. | | |
| Qt./Gal.†..... | Conductivity/TDS Standard Solution (1 ltr./3,8 ltr.)..... | \$29.50/103.50 |
| Gallon†..... | Custom Conductivity/TDS Standard Solution (3,8 ltr.)..... | NET \$162.00 |
| 2 oz.†..... | Conductivity/TDS Standard Solution (59 ml) (Available for 442-3000, 442-1500, KCl-7000, KCl-1800, NaCl-14.0 and NaCl-7500 only)..... | \$16.00 |
| pH Qt./Gal.†..... | pH Calibration Buffer Solution 4, 7 or 10 (1 ltr./3,8 ltr.)..... | \$29.50/103.50 |
| pH 2 oz.†..... | pH Calibration Buffer Solution 4, 7 or 10 (59 ml)..... | \$16.00 |
| SSQ/SSG..... | pH/ORP Sensor Storage Solution (1 qt./ltr. or 1 gal./3,8 ltr.)..... | \$29.50/103.50 |
| SS2OZ..... | pH/ORP Sensor Storage Solution (2 oz./59 ml)..... | \$16.00 |
| † IMPORTANT: Exact value and size (e.g., 442-3000Q or PH7G) must be specified when ordering solutions. Above prices do not include a written certificate; see item SC below. | | |

NIST CERTIFICATIONS - NIST (National Institute of Standards & Technology) Certification

All conductivity and pH solutions used in factory calibration are NIST traceable. Written certification of factory NIST calibration of instruments and/or NIST traceability of conductivity and pH solutions are available at the prices below. **Note:** Test data not included.

Certificate of factory NIST calibration with data for DIGITAL DIALYSATE METER

| | | |
|------------|--|--------------|
| MCD..... | Ordered <i>with</i> meter - Must specify instrument model, e.g., MCD-D-6..... | NET \$144.00 |
| SC-XX..... | Ordered <i>with</i> solution - Must specify solution values, e.g., SC-3000Q..... | NET \$57.00 |
| SCATF..... | Ordered <i>after</i> solution is shipped*..... | NET \$96.50 |

* Serial number, or batch number and expiration date must be supplied when ordering.

TERMS AND CONDITIONS (Additional information upon request)

All special order items are non-returnable and non-refundable.

Minimum order \$25.00. All prices are US dollars and are F.O.B. Carlsbad, CA.

Terms: 1% 10 days, net 30 days on approved credit.

Credit Cards accepted: Visa, MasterCard, American Express, and Discover.

Because of our policy of continuous product improvement, the Myron L® Company reserves the right to make changes in design, specifications, and prices without notice.

2450 Impala Drive
Carlsbad, CA 92010-7226 USA
Tel: +1-760-438-2021
Fax: +1-800-869-7668/+1-760-931-9189
www.myronl.com

Built On Trust. Founded in 1957, the Myron L® Company is one of the world's leading manufacturers of water quality instruments. Because of our commitment to product improvement, changes in design and specifications are possible. You have our assurance any changes will be guided by our product philosophy: accuracy, reliability, and simplicity.

MYRON L®
COMPANY
Water Quality Instrumentation
Accuracy • Reliability • Simplicity